



Water Bath Shaker

LB-66WBS

Index

Sr.no	Title	Page no
1.	Safety Measures	2
2.	Introduction	3
3.	Features	3
4.	Specifications	4
5.	Applications	4
6.	Operations	5

1. Safety Measures

- **Timing Adjustment:** Set the timer as needed; if not required, turn it to "ON."
- **Handling Vessels:** Always turn off the instrument before placing or removing a vessel.
- **Stable Placement:** Ensure the instrument is on a stable, level surface.
- **Grounding & Power Off:** Use a properly grounded outlet and turn off the power after use.
- **Product Maintenance:** Our products are continuously maintained. Improper or intentional misuse may incur additional costs.

2. Introduction

Water Bath Shaker LB-66WBS has a 31-liter capacity and operates across a temperature range from RT to 100°C for accurate heating. It offers a rotation speed range of zero to 300 rpm with precision of minus plus ± 1 rpm ensuring variable mixing. It accommodates a range of flask size, providing flexibility for various sample volumes. Our Shaker is constructed with durable stainless steel ensuring the reliability across various scientific field.

3. Features

- 1) Microprocessor control unit
- 2) Stainless steel spring clamp for multipurpose use
- 3) Constructed from high quality material
- 4) Even heat distribution with water circulation
- 5) Offer over temperature protection

4. Specifications

Model	LB-66WBS
Capacity	31 L
Temperature range	RT to 100
Rotation speed range	0 to 300 rpm
Speed position	± 1 rpm
Timer range	0 to 120 min
Flask Capacity	24×100 ml 12×250 ml 6×500 ml 4×1000 ml
Oscillation modes	3 votes rotation reciprocating and continuous
Display	20 mm
Power consumption	LED
Power supply	1550 W
Chamber dimension	AC 110/220 V, 50/60 Hz
Product dimension	390 × 490 × 165 mm
Packaging dimension	720 × 460 × 520 mm
Lightweight net weight	50 kg
Gross weight	61 kg

5. Applications

What about shakeup LB-66WBS is used for temperature controls shaking in laboratory experiments research application and pharmaceutical process

6. Operations

1) Power Supply Requirements:

Ensure that the instrument is connected to an AC 220V power source. The outlet must have a current rating of at least 6A and be properly grounded for safety. If the outlet does not meet these requirements, replace it before use.

2) Setting Up the Instrument:

- Open the instrument cover.
- Secure the vessel requiring oscillation within the spring multipurpose clamp.

3) Operation:

- Turn on the power switch and the oscillation switch.
- Select the desired oscillation speed and operating temperature.
- Gradually adjust the oscillation speed from slow to fast as needed.
- Set the required temperature using the temperature control unit.

4) Indicator Lights:

- **Green Light:** Indicates that heating is active.
- **Red Light:** Indicates that heating is stopped.
- The displayed temperature represents the actual working chamber temperature.



Labotronics Scientific. 1007 N Orange St., 1382, Wilmington, DE 19801, USA.
info@labotronics.com | www.labotronics.com